



淮安超越医疗用品厂

CHAOYUE MEDICAL INSTRUMENTS FACTORY

CHAO YUE MEDICAL manufactures high quality Single-Use stainless steel ophthalmic micro-surgical instruments , plastic surgical instruments , Neuro-surgical instruments in its modern factory complex in Huai'an, China . With more than 15 years experience, we have upgraded our production facility from handle-making to semi-auto production which have greatly improved the quality consistency and delivery .

Quality and customer satisfaction are the foremost, we have the compressive quality management procedure with qualified staffs and high-tech equipments. we have also been certified with ISO13485 by NSF-ISR .

We will be your trusted partner in single-use stainless steel surgical instruments .



TEL: +86 18970878605 13852367707

E-mail: sales@chaoyue-medical.com



| | |
|----------------------------------|-------|
| 1 Forceps | 1-13 |
| 2 Speculum | 14-16 |
| 3 Needle-Holders | 17-18 |
| 4 Sprung Scissors | 19-20 |
| 5 Clamps & Curettes | 21-22 |
| 6 Indentors & Depressors | 23 |
| 7 Markers & Gauges | 24-25 |
| 8 Lacrimal | 26-27 |
| 9 LASIK | 28-29 |
| 10 Choppers & Manipulators | 30-37 |
| 11 Miscellaneous | 38 |

1 Capsulorhexis Forceps

Cross Action




Short Handled Direct Action




Long Handled Direct Action


Sharp Tip




Blunt Tip




Curved Direct




Straight Direct



Curved





Straight



Fixation Forceps 2



CYSU1114

Fine Notched Forceps, Short
(Hoskins)
Ø0.2mm Notch Tip,
5mm Tying Platforms.
52mm Handle




CYSU1115

Notched Forceps
(Hoskins)
Ø0.3mm Notch Tip,
5mm Tying Platforms.
80mm Handle




CYSU1123

Fine Toothed Forceps, Short
(Bonn)
0.3mm 1×2 Teeth,
5mm Tying Platforms.
52mm Handle



CYSU1124

Toothed Forceps (Bonn)
0.4mm 1×2 Teeth,
5mm Tying Platforms.
80mm Handle




Long-Handled

Short-Handle




CYSU1209

Colibri Notched Forceps, Short
Ø0.2mm Notch Tip,
5mm Tying Platforms.
52mm Handle




CYSU1128

Colibri Toothed Forceps
0.4mm 1×2 Teeth,
5mm Tying Platforms.
72mm Handle




CYSU1129

Colibri Toothed Forceps, Short
0.4mm 1×2 Teeth,
5mm Tying Platforms.
52mm Handle




CYSU1168

Toothed Forceps, 15mm Shaft
0.4mm 1×2 Teeth,
15mm Shaft
7.2mm Tying Platforms.
80mm Handle




3 Capsulorhexis Forceps

Direct Action MICS Capsulorhexis Forceps




CYSU1567

MICS Capsulorhexis - Curved/Sharp
 0.4mm × 0.3mm Sharp Tips,
 12.0mm Curved Direct Action Jaws,
 1.2mm at base of Shaft.
 80mm Flat Handle




CYSU1568

MICS Capsulorhexis - Straight/Sharp
 0.4mm × 0.3mm Sharp Tips,
 12.0mm Straight Direct Action Jaws,
 1.2mm at base of Shaft.
 80mm Flat Handle




Flat Handled MICS Cross-Action Capsulorhexis Forceps




CYSU1236

MICS Capsulorhexis - Curved/Sharp
 0.4mm × 0.3mm Sharp Tips,
 14.5mm Curved Cross-Action Jaws,
 1.0mm × 1.65mm at Hinge.
 27.5mm Flat Handle



CYSU1237


MICS Capsulorhexis - Straight/Sharp
 0.4mm × 0.3mm Sharp Tips,
 14.5mm Straight Cross-Action Jaws,
 1.0mm × 1.65mm at Hinge.
 27.5mm Flat Handle



Capsulorhexis Forceps 4


CYSU1121

Capsulorhexis - Curved/Sharp
 0.6mm × 0.3mm Sharp Tips,
 12mm Curved Direct Action Jaws.
 80mm Flat Handle




CYSU1122

Capsulorhexis - Straight/Sharp
 0.6mm × 0.3mm Sharp Tips,
 12mm Straight Direct Action Jaws.
 80mm Flat Handle




CYSU1164

Capsulorhexis - Curved
 0.4mm × 0.3mm Blunt Tips,
 12mm Curved Direct Action Jaws.
 80mm Flat Handle




CYSU1165

Capsulorhexis - Straight
 0.4mm × 0.3mm Blunt Tips,
 12mm Straight Direct Action Jaws.
 80mm Flat Handle




Long Handled Direct Action Capsulorhexis Forceps




CYSU1166

Capsulorhexis - Short/Curved/Sharp
 0.4mm × 0.3mm Sharp Tips,
 12mm Curved Direct Action Jaws.
 53mm Flat Handle




CYSU1167

Capsulorhexis - Short/Straight/Sharp
 0.4mm × 0.3mm Sharp Tips,
 12mm Straight Direct Action Jaws.
 53mm Flat Handle




Short Handled Direct Action Capsulorhexis Forceps



5 Capsulorhexis Forceps


CYSU1135

Capsulorhexis - X-Action/Curved
 0.4mm × 0.3mm Blunt Tips,
 12mm Curved Cross-Action Jaws.
 27.5mm Flat Handle




CYSU1234

Capsulorhexis - X-Action/Straight
 0.4mm × 0.3mm Blunt Tips,
 12mm Straight Cross-Action Jaws.
 27.5mm Flat Handle




CYSU1136

Capsulorhexis - X-Action/Curved/Sharp
 0.4mm × 0.3mm Sharp Tips,
 12mm Curved Cross-Action Jaws.
 27.5mm Flat Handle




CYSU1235

Capsulorhexis - X-Action/Straight/Sharp
 0.4mm × 0.3mm Sharp Tips,
 12mm Straight Cross-Action Jaws.
 27.5mm Flat Handle




Flat Handled Cross Action Capsulorhexis Forceps




CYSU1240

MICS Capsulorhexis - Round/Curved/Sharp
 0.4mm × 0.3mm Sharp Tips,
 11mm Curved Cross-Action Jaws,
 1.0mm × 1.65mm at Hinge.
 29mm Round Handle




CYSU1241

MICS Capsulorhexis - Round/Straight/Sharp
 0.4mm × 0.3mm Sharp Tips,
 11mm Straight Cross-Action Jaws,
 1.0mm × 1.65mm at Hinge.
 29mm Round Handle




Round Handled MICS Cross-Action Capsulorhexis Forceps



Fixation Forceps 6


CYSU1139

45° Toothed Forceps
0.3mm 1×2 Teeth,
5mm Tying Platforms.
80mm Handle



CYSU1145

45° Fine Toothed Forceps, Short
0.27mm 1×2 45° Teeth,
5mm Tying Platforms.
52mm Handle



CYSU3010

Converse Green Fixation Forceps
2 × 7.9mm Tips with 0.1mm × 1.0mm Teeth.
102mm Handle




CYSU1442

Corneal Fixation Forceps
2 Pairs of 0.15mm × 0.16mm 1×2 Teeth,
12.2mm Apart on Curved Shaft,
Angled at 50°.
52mm Handle



CYSU1495

Fixation Ring
0.5mm Thick 13mm Internal Diameter Ring,
1.5mm × Ø0.6mm Teeth at 3.5mm intervals,
1.3mm × Ø0.4mm Teeth at 0.9mm intervals.
75mm Handle



CYSU1137

Thornton Fixation Ring
0.5mm Thick 6.5mm Radius Ring
2mm x Ø0.9mm Teeth at 5mm intervals.
75mm Handle



7 Tying Forceps

CYSU1125

Straight Tying Forceps
0.35mm Tip,
5mm Tying Platforms.
80mm Handle



CYSU1219

Straight Tying Forceps, Short Handle
0.26mm Tip,
5mm Tying Platforms.
52mm Handle



CYSU1210

Curved Tying Forceps
0.38mm Tips,
10mm Curved Tying Platforms.
80mm Handle



CYSU1127

Angled Tying Forceps
0.35mm Tip,
5mm Tying Platforms,
8mm Shaft Angled 45°.
80mm Handle



CYSU1413

Catford Forceps
0.55mm × 0.8mm Tips,
7.0mm Triangular Shaped Serrated Platforms.
73mm Long



CYSU1414


Birks Tying Forceps
0.1mm × 0.3mm Tips,
5mm Tying Platforms.
77mm Round Handle



Lens Forceps 8


CYSU1126

IOL Implantation Forceps
7mm Highly Polished Curved Convex Jaws.
80mm Handle



CYSU1138

Lens Folding Forceps
Rounded Inner Jaws,
9.5mm 45° Angled Shaft.
80mm Handle




CYSU1131

Lens Folding Forceps
Triangular Shaped Inner Jaws,
10.5mm 45° Angled Shaft.
80mm Handle



CYSU1460

IOL Removal Forceps
0.1mm × 0.1mm Ridges Along 4mm of 15mm
Shaft, Angled at 30°.
52mm Handle



CYSU1146

Lens Loading Forceps
0.6mm × 8mm Tips,
45° Angled Shaft.
80mm Handle



CYSU1154

ICL Lens Loader Forceps
0.6mm × 8mm Tips,
45° Angled Shaft.
80mm Handle



CYSU1284

Contact Lens Forceps
0.6mm × 0.4mm Tips,
9mm Curved Platforms.
70mm Handle



CYSU1169

LASIK Flap-Lifting Forceps
0.3mm × 11.5mm Tips,
45° Angled Shaft,
80mm Handle.



9 General Forceps

CYSU0277

Moorfields Forceps
1mm × 1.25mm Tips,
13.5mm Serrated Jaws.
70mm Handle




CYSU0247

Iris Pointed Forceps, 118mm
0.5mm × 0.6mm Tips,
14mm Serrated Jaws.
70mm Handle




CYSU0278

Iris Fine Toothed Forceps
1.2mm × 2.0mm 1×2 Teeth.
70mm Handle




CYSU0080

Fine Stem Toothed Forceps
2mm × 1mm 1×2 Teeth.
70mm Handle



CYSU0274

Jayles Forceps
0.3mm 1×2 Teeth
Angled Forward 40°.
70mm Handle




CYSU1445

Maxfine Forceps
0.2mm x 0.2mm x 11mm Curved Tips.
70mm Handle




CYSU0276

Bishops Forceps
0.7mm × 0.9mm 1×2 Teeth.
70mm Handle



CYSU0279

St Martins Forceps
0.90mm × 0.26mm 1×2 Teeth.
53mm Handle



CYSU0085

Epilation Forceps
5mm × 7mm Flat Plane Tips.
70mm Handle













CYSU1418

Meibomian Gland Forceps
3.0mm × 3.5mm x 1.0mm Block Tips
3.5mm 35° Angled Shaft.
70mm Handle



Jewellers Forceps 10

| | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <p>CYSU0286 Jewellers Forceps, Fine Straight 7mm Pointed Tip. 80mm Handle</p> |  <p>CYSU0283 Jewellers Forceps, Fine Angled 7mm 45° Angled Pointed Tip. 80mm Handle</p> |
|  <p>CYSU1472 Jewellers Forceps, Angled/Smooth 7mm 45° Angled Pointed Tip. 80mm Handle</p> |  <p>CYSU0282 Jewellers Forceps, Curved/Serrated 0.5mm Fine Tip Curved and Serrated 15mm Jaws. 80mm Handle</p> |
|  <p>CYSU0284 Jewellers Forceps, Curved/Plain Fine Point, "Swan Neck" Curved Shaft. 85mm Handle</p> |  <p>CYSU0238 Adson Forceps, Toothed 1.5mm × 2.75mm 1×2 Teeth. 46mm Handle</p> |
|  <p>CYSU0241 Adson Forceps, Non-Toothed 1.5mm × 12.5mm Serrated Jaws. 46mm Handle</p> |  <p>CYSU0075 Punctum Plug Forceps 1.8mm x 2.6mm 45 Degree Angled Tips with 0.8mm slot. 84mm Handle</p> |
|  <p>CYSU1411 Scleral Plug Forceps 1.15mm × 2.3mm Tips with Ø0.5mm Notch, 6mm 55° Angled Shaft. 65mm Handle</p> |  <p>CYSU1498 Pre-Chopper Forceps 4.5mm × 1.6mm Tips with 3.0mm Spread, 34.5mm Straight Cross-Action Jaws. 73mm Flat Handle</p> |

11 Specula - Mechanical



CYSU1258

Clarke Speculum
8.5mm × 14.5mm Open Blades,
Ø1.5mm Wire.
Metal Body and Screw
Mechanical Speculum




CYSU1263

Bangerter Speculum, Open Blade
5.0mm × 10.5mm Open Blades,
Ø1.0mm Wire.
Metal Body and Screw
Mechanical Speculum



CYSU1267

Bangerter Speculum, Solid Blade
6.0mm × 13.5mm Solid Blades,
0.7mm Thick.
Metal Body and Screw
Mechanical Speculum



CYSU1266

Barrett Speculum
6mm × 14.5mm Curved Blades,
32mm Wide
Screw Locking Mechanism

Specula - Sprung 12


CYSU1287

Kratz-Barraquer Speculum, Curved
6.5mm × 14.5mm Open Curved Blades,
Ø0.9mm Wire.
30mm Wide




CYSU1252

Kratz-Barraquer Speculum
6.5mm × 14.0mm Open Blades,
Ø0.9mm Wire.
30mm Wide



CYSU1270

Pierce Adjustable Speculum
6.5mm × 12.5mm Open Blades,
Ø0.8mm Wire.
29mm Wide Adjustable
to 48.5mm Wide



CYSU1275

Barraquer Speculum, Flat Blade
6.5mm × 14.5mm Closed Blades,
Ø1mm Wire.
30mm Wide



CYSU1250

Barraquer Speculum, Curved Blade
5.0mm × 16.5mm Curved Closed Blades,
Ø1.0mm Wire.
38mm Wide













CYSU1264

Adjustable Speculum
6.5mm × 14.0mm Open Blades,
Ø0.9mm Wire.
25mm Wide Adjustable
to 45mm Wide



13 Specula-Neo Natal/Paediatric

| | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>CYSU1265 Paediatric Mechanical Speculum, Open 5.0mm × 10.5mm Open Blades, Ø1.0mm Wire. Polycarbonate Body and Screw Mechanical Speculum</p>  | <p>CYSU1560 Mechanical Paediatric Speculum, Closed 5.0mm × 10.5mm Closed Blades, Ø1.0mm Wire. Polycarbonate Body and Screw Mechanical Speculum</p>  |
| <p>CYSU1260 Alphonso Paediatric Speculum 7.0mm × 5.5mm/8.0mm × 8.5mm Closed Blades, Ø1.0mm Wire. 48mm Wide</p>  | <p>CYSU1261 Neonatal Speculum, Large 4.0mm × 10mm Closed Blades, Ø0.8mm Wire. 25mm Wide</p>  <p>CYSU1262 -Neonatal Speculum, Small 3.0mm × 7.0mm Closed Blades</p> |
| <p>CYSU1278 Paediatric Speculum, 3cm 5.5mm × 13.5mm Closed Blades, Ø0.9mm Wire. 31.5mm Wide</p>  | <p>CYSU1292 Neonatal Speculum, Solid Blade 5.0mm × 5.0mm Solid Blades, 0.1mm Thick 22mm Wide</p>  |
| <p>CYSU1273 Cook Neonatal Speculum 5.5mm × 8.0mm Solid Curved Blades, 0.1mm Thick, Locking Screw. Adjustable to 38mm Wide</p>  | <p>CYSU1259 Sauer Neonatal Speculum 4.6mm × 6.6mm Solid Blades, Ø0.6mm Wire Spring Arms. 18.5mm Wide</p>  |
| <p>CYSU1293 Infant Speculum 4.3mm × 4.3mm Closed Blades, Ø0.7mm Wire. 25mm Wide</p>  | <p>CYSU1288 Corcelle Speculum, Solid Blade 6.5mm × 8.5mm Solid Blades, 0.7mm Thick. 55mm Wide</p>  |

Needle-Holders 14



CYSU1110

Corneal Needle-Holder
4mm Fine Curved Jaws,
0.45mm Tip.
32mm Handle



CYSU1111

Corneal Needle-Holder, Round Handle
4mm Fine Curved Jaws,
0.45mm Tip.
118mm Round Handle



CYSU3015

Barraquer Needle-Holder
0.42mm × 0.45mm Tips,
4mm Fine Curved Jaws with 5mm Shaft.
138mm Long



CYSU3017

Corneal Needle-Holder, Long
0.45mm Tip, 9mm Curved Jaws,
Round Handle.
138mm Long

15 Needle-Holders



CYSU1243

Troutman Needle-Holder
4mm Fine Curved Jaws,
0.55mm Tip.
112mm Round Handle



CYSU3009

Castroviejo Needle-Holder, Without Lock
1.5mm Tips, 12.5mm Straight Jaws,
Non-Locking.
142mm Long



CYSU3011

Castroviejo Needle-Holder, Without Lock
1.5mm Tips, 12.5mm Curved Jaws,
Non-Locking.
142mm Long



CYSU3012

Castroviejo Needle-Holder, Straight
1.5mm Tips,
12.5mm Serrated Cross-Action Jaws With Lock.
142mm Long



CYSU3013

Castroviejo Needle-Holder, Curved
1.5mm Tips,
12.5mm Serrated Cross-Action Jaws With Lock.
142mm Long

Sprung Scissors 16



CYSU1211

Westcott Scissors, Curved, Pointed
0.5mm × 0.6mm Pointed Tips,
17.5mm Curved Blades.
36mm Handle



CYSU1212

Westcott Scissors, Straight, Rounded
0.5mm × 0.6mm Rounded Tips,
17.5mm Straight Blades.
36mm Handle



CYSU1213

Westcott Scissors, Curved, Dolphin Nose
17.5mm Curved 'Dolphin Nose' Blades.
36mm Handle



CYSU1214

Westcott Scissors, Straight, Pointed
Pointed Tips,
17.5mm Straight Blades.
36mm Handle



CYSU1217

Castroviejo Scissors, Straight, Rounded
0.5mm × 0.6mm Rounded Tips,
11.5mm Straight Blades.
36mm Handle




CYSU1218

Castroviejo Scissors, Curved, Rounded
0.5mm × 0.6mm Rounded Tips,
11.5mm Curved Blades.
34mm Handle



17 Sprung Scissors



| | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <p>CYSU1140 Vannas Scissors, Curved Pointed Tips on 8mm Curved Blades, Ø1.5mm When Closed. 34mm Handle</p> |  <p>CYSU1141 Vannas Scissors, Straight Pointed Tips on 8mm Straight Blades, Ø1.5mm When Closed. 34mm Handle</p> |
|  <p>CYSU1469 Castroviejo Corneoscleral Scissors (L) 0.35mm x 0.3mm Pointed Tips, 7mm Curved Blades. 34mm Handle</p> |  <p>CYSU1470 Castroviejo Corneoscleral Scissors (R) 0.35mm x 0.3mm Pointed Tips, 7mm Curved Blades. 34mm Handle</p> |
|  <p>CYSU1116 Ong Capsule Scissors Pointed Tips, 10mm Angled Blades, 2.2mm x 2.3mm Thick at Hinge. 34mm Handle</p> |  <p>CYSU1459 IOL Cutter (Osher) 0.2mm x 0.5mm Blunt Tips 9.5mm Serrated Blades. 66mm Handle</p> |

Clamps & Curettes 18



CYSU1148

Chalazion Clamp, Oval, Large
15mm × 25mm Oval Profile,
12mm × 22mm Aperture.
90mm Long

CYSU1149

Chalazion Clamp, Oval, Medium
13mm × 20mm Oval Profile,
10mm × 17mm Aperture.
90mm Long

CYSU1150

Chalazion Clamp, Round
Ø12.0mm Round Profile,
Ø9mm Round Aperture.
90mm Long

CYSU1448

Chalazion Clamp, Oval, Large, Grooved
15mm × 25mm Oval Profile,
Grooved, 12mm × 22mm Aperture.
90mm Long

CYSU1449

Chalazion Clamp, Oval, Medium, Grooved
13mm × 20mm Oval Profile, Grooved,
10mm × 17mm Aperture.
90mm Long

CYSU1450

Chalazion Clamp, Round, Grooved
12.0mm Round Profile, Grooved,
Ø9mm Round Aperture.
90mm Long

19 Clamps & Curettes



| | | | |
|--|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| | <p>CYSU1157 Chalazion Curette, Double, Ø1mm/Ø2mm Ø1.0mm Curette, 1.0mm Depth Ø2.0mm Curette, 1.35mm Depth. 128mm Long</p> | <p>CYSU1158 Chalazion Curette, Double, Ø2mm/Ø3mm Ø2.0mm Curette, 1.35mm Depth Ø3.0mm Curette, 1.6mm Depth. 128mm Long</p> | |
|--|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|

CYSU1151
 Chalazion Curette, Ø1mm
 Ø1.0mm Curette,
 1.0mm Depth.
 115mm Long

CYSU1155
 Chalazion Curette, Ø2mm
 Ø2.0mm Curette,
 1.35mm Depth.
 115mm Long

CYSU1156
 Chalazion Curette, Ø3mm
 Ø3.0mm Curette,
 1.6mm Depth.
 115mm Long

CYSU1159
 Chalazion Curette, Ø4mm
 Ø4.0mm Curette,
 1.7mm Depth.
 115mm Long



Indentors & Depressors 20

CYSU1405

Schocket Double-Ended Scleral Depressor

2.7mm × 6.5mm T Shaped Tip,

5.5mm × 2.2mm Rounded Tip.

68mm Handle



CYSU1407

VR Double-Ended Depressor/Indentor

2.7mm × 6.5mm T Shaped Tip,

5.5mm × 2.2mm Rounded Tip.

68mm Handle



CYSU1409

Curved Scleral Depressor

Ø1.5mm × 6.1mm T Shaped Tip,

21mm 50° Angled Curved Shaft.

78mm Handle



CYSU1247


Triangular Scleral Depressor

7.0mm × 7.0mm Triangular Shaped Tip,

Ø0.75mm,

15mm Straight Shaft.

76mm Handle



CYSU1408

Ball Tip Scleral Depressor

Ø3.0mm Ball Tip,

21.5mm Straight Shaft.

50mm Handle



21 Markers & Gauges



CYSU1401

Castroviejo Caliper, Straight Tip

0.4mm × 0.8mm Straight Tips,

30mm Long.

0mm to 25mm in 0.5mm Gradations




CYSU1403

Castroviejo Caliper, Angled Tip

0.25mm × 0.3mm Angled Tips,

34mm Long.

0mm to 25mm in 0.5mm Gradations




CYSU1229

Braunstein IVT Marker, 3.5mm/4.0mm

2 × 1.5mm Tips with 3.5mm Spread,

2 × 1.5mm Tips with 4.0mm Spread.

90mm Long



Markers & Gauges 22


CYSU1195

Pointed IVT Marker, 3.5mm/4.0mm

2 × Ø0.95mm Tips with 3.5mm Spread,

2 × Ø0.95mm Tips with 4.0mm Spread.

108mm Long




CYSU1286

Pointed IVT Marker 2.5/3.0mm

2 × Ø0.95mm Tips with 2.5mm Spread,

2 × Ø0.95mm Tips with 3.0mm Spread.

108mm Long



CYSU1285

Double Ended IVT Guide Pad, 3.5/4.0mm

Ø7.0mm Pad with 3.5mm Guide, 6.0mm 45° Angled Shaft,

Ø0.8mm Pad with 4.0mm Guide, 6.0mm 45° Angled Shaft.

100mm Long



23 Lacrimal



CYSU1417

Golf Club Spud

0.1mm x 0.8mm x 2.0mm Oval Profile 'Golf Club',

18mm Straight Shaft.

75mm Handle



CYSU1180

Punctum Dilator, Medium

Ø0.55mm x 24mm Tip,

76mm Handle.



CYSU1181

Punctum Dilator, Fine/Medium

Ø0.4mm x 21.5mm Tip

Ø0.6mm x 18mm Tip,

77mm Handle.



CYSU1183

Punctum Seeker

5mm Pointed Tip, 18mm Straight Shaft,

77mm Handle.



CYSU1193

Punctum Dilator/Seeker

Ø0.45mm x 21.5mm Tip

5mm Pointed Tip, 18mm Straight Shaft.

77mm Handle.



Lacrimal 24


CYSU1426

Punctum Dilator, Small
Ø0.3mm × 20mm Dilator.
30mm Handle.



CYSU1427

Lacrimal Probe
1.0mm Pointed Tip, Ø0.6mm 80mm Probe.
15mm Handle.



CYSU1419

Punctum Dilator/Lacrimal Probe
Ø0.28mm × 25mm Dilator
Ø0.7mm × 44.5mm Probe.
70mm Handle.



CYSU1420

Bowmans Probes, Size 00-000
Ø0.8mm (00g) × 50mm Probe
Ø0.7mm (000g) × 50mm Probe.
30mm Handle.




CYSU1421

Bowmans Probes, Size 0-00
Ø0.9mm (0g) × 50mm Probe
Ø0.8mm (00g) × 50mm Probe.




















CYSU1422

Bowmans Probes, Size 1-2
Ø1.0mm (1g) × 50mm Probe
Ø1.1mm (2g) × 50mm Probe.



25 LASIK

| | | | |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  | <p>CYSU1177 LASIK Flap Elevator, Curved, 8mm 0.2mm × 0.15mm Tip, Curved 8mm Tapered Shaft at 90°. 75mm Handle</p> | <p>CYSU1178 LASIK Flap Elevator, Straight, 8mm 0.15mm × 0.4mm Flattened Tip, 8mm 45° Angled Shaft. 75mm Handle</p> |  |
|  | <p>CYSU1435 LASIK Flap Elevator, Curved, 10mm 0.2mm × 0.15mm Tip, Curved 10mm Tapered Shaft at 90°. 75mm Handle</p> | <p>CYSU1206 LASIK Flap Separator & Elevator (Binder) 0.1mm × 0.3mm Rounded Tip 0.7mm Long, 10mm Curved 60° Angled Shaft. 75mm Handle</p> |  |
|  | <p>CYSU1133 LASIK Flap Separator/Elevator 1mm Micro-Finger Tip, 11mm Curved Shaft, Angled 45°. 75mm Handle</p> | <p>CYSU1171 LASIK Flap Elevator, Pointed Hook Ø0.5mm × 11.5mm Pointed Hook. 75mm Handle</p> |  |
|  | <p>CYSU1175 LASIK Spatula (Johnson) 0.5mm × 1.0mm Pointed Tip, 12mm Curved 45° Angled Shaft. 75mm Handle</p> | <p>CYSU1184 LASIK Flap Elevator, Beam 0.35mm × 10mm Beam Between 18mm Shafts. 75mm Handle</p> |  |
|  | <p>CYSU1176 LASIK Flap Separator, Flattened Tip 0.75mm × 1.1mm Flattened Tip, 23mm Straight Shaft. 75mm Handle</p> | <p>CYSU1189 Double Ring Corneal Marker r1.75mm Ring Overlapping r2.0mm Ring, 2mm Deep with Ridged Underside. 75mm Handle</p> |  |
|  | <p>CYSU1134 "Y" Shape Corneal Marker 2 × 8mm Marking Blades. 75mm Handle</p> | | |

| | |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <p>CYSU1172 LASIK Flap Separator/Elevator (Seibel) 0.16mm × 0.28mm × 1.0mm Tip, 20mm Straight Shaft. Ø0.5mm × 11.5mm Pointed Hook. 52mm Handle</p> |  <p>CYSU1188 LASIK Flap Elevator/Shield 0.1mm × 1.0mm Tip, 11mm 45° Curved Angled Shaft. 0.72mm Lip, 3.5mm × 8.0mm Shield, 3.5mm 45° Angled Shaft. 52mm Handle</p> |
|  <p>CYSU1205 LASIK Flap Separator/Elevator 2.95mm × 8mm Triangular Tip at a 60° Angle, 24mm Straight Shaft. Ø0.53mm Pointed Tip, 9mm 50° Angled Shaft. 52mm Handle</p> |  <p>CYSU1227 LASIK Cap/Flap Elevator 0.4mm × 1.2mm Rounded Tip, Curved Shaft 0.4mm × 1.8mm Flat Tip, 0.65mm Curved Shaft. 52mm Handle</p> |
|  <p>CYSU1294 LASIK Flap Separator/Elevator, Curved 0.15mm × 0.30mm × 1.0mm Separator. 0.15mm × 0.30mm × 8.0mm Curved Elevator at 90° to Shaft. 52mm Handle</p> |  <p>CYSU1295 LASIK Flap Separator/Elevator, Straight 0.15mm × 0.30mm × 1.0mm Separator. 0.15mm × 0.30mm × 8.0mm Straight Elevator at 90° to Shaft. 52mm Handle</p> |

CYSU1296

Paton Spatula/Spoon

2.0mm × 7.5mm Spatula Rounded Tip, 0.25mm Thick at 30° Angle to 24.0mm Shaft

5.0mm × 7.0mm Spoon Tip, 0.3mm Thick, 26.0mm Shaft.

52mm Handle




27 Choppers & Manipulators

CYSU1437



CYSU1238



CYSU1142



CYSU1179



CYSU1108



Hooks & Retractors 28













| | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <p>CYSU3005 Guthrie Hook Delicate Hooks 2.0mm Wide</p> |  <p>CYSU3006 Gillies Skin Hook, 18cm 2.0mm × 2.5mm Hook, 77mm Shaft. 180mm Long</p> |
|  <p>CYSU10232 Chevasse Muscle Hook, 'W' Shaped 12mm 'W' Shaped Hook at 90°. 75mm Handle</p> |  <p>CYSU0253 Curved Muscle Hook 10.3mm Curved Hook. 75mm Handle</p> |
|  <p>CYSU3001 Desmarres Lid Retractor, Small 11.0mm × 11.0mm × 5.0mm Retractor Tip, 35mm Straight Shaft. 76mm Handle</p> |  <p>CYSU0252 Ball Tip Muscle Hook 10mm Ball Tip Hook. 90mm Handle</p> |
|  <p>CYSU3038 Fisons Retractor 10.8mm × 20.0mm Flat Retractor Blade, 3.0mm Lip. 76mm Handle</p> |  <p>CYSU1444 Lens Expression Loop Ø0.4mm Wire with 0.1mm Ridges Forming 1.5mm × 13mm Loop. 76mm Handle</p> |

CYSU1425













Eye Guard
18.0mm × 6.0mm Rounded End
20.5mm × 6.5mm Rounded End.
100mm Long















29 Choppers & Manipulators

| | | | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <p>CYSU1480</p> <p>Reverse Dialling Hook 0.92mm Cylindrical Tip, Rounded Bottom 9.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1489</p> <p>Neuhann Nucleus Divider 0.5mm Angled Rounded Tip, 9.5mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1105</p> <p>Claw Chopper 0.3mm × 1.5mm Hooked Tip, Rounded Bottom, 11.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1446</p> <p>Claw Chopper, Fine 0.35mm × 1.5mm Hooked Tip, Rounded Bottom, 12.0mm 40° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1104</p> <p>Sadiq Chopper 0.4mm × 1.7mm Curved Tip, Rounded Bottom, 10mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1106</p> <p>Rosen Chopper 0.35mm × 1.2mm Straight Tip, Rounded Back, 10mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1119</p> <p>Greens Chopper, Right Handed 0.3mm × 1.4mm Chopper Tip, Rounded Bottom, 10mm 45° Angled Shaft, 20° Bend. 75mm Handle</p> |  | <p>CYSU1120</p> <p>Greens Chopper, Left Handed 0.3mm × 1.4mm Chopper Tip, Rounded Bottom, 10mm 45° Angled Shaft, 20° Bend. 75mm Handle</p> |  |
| <p>CYSU1108</p> <p>Phaco Chopper 0.3mm × 1.5mm Chopper Tip, Rounded Bottom, 11.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1143</p> <p>McPherson Manipulator 0.15mm × 0.4mm × 1.7mm Flattened Tip , Angled at 30° from 10.6mm, 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1107</p> <p>Jacobs Chopper 2mm Hook Tip, 20mm Straight Shaft. 75mm Handle</p> |  | <p>CYSU1569</p> <p>Femto Manipulator 0.2mm × 0.7mm × 13mm 60° Angled Shaft. 75mm Handle</p> |  |

Choppers & Manipulators 30

| | | | |
|-----------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <p>CYSU1144 Drysdale Manipulator 0.95mm × 1.6mm 'Rudder' Tip, 10mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1226 Spalton Chopper 1.7mm Hooked Chopper Tip, 11mm Curved 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1441 Nagahara Chopper 0.3mm × 1.4mm Wedge Tip, 10mm 50° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1118 Koch Chopper, 1.9mm 1.9mm Cylindrical Tip, Rounded Bottom, 9.5mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1228 Koch Chopper, 1.4mm 1.4mm Cylindrical Tip, Rounded Bottom, 9.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1440 Koch Chopper, 1.4mm, Thin Shaft 1.4mm Cylindrical Tip, Rounded Bottom, Ø0.25mm × 9.5mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1471 Denis Nucleus Chopper 0.3mm Curved Cleaver Tip Angled Left 30° from 11mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1103 Sinskey Hook, Ø0.2mm Ø0.2mm × 0.5mm Hook, 10mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1173 Sinskey Hook, Ø0.1mm Ø0.1mm × 0.6mm Hook, 10mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1147 Sinskey Hook, Straight Ø0.18mm × 0.5mm Hook, 16mm Straight Shaft. 75mm Handle</p> |  |
| <p>CYSU1436 Hagège Manipulator Ø0.2mm × 0.5mm Hook, 11.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1102 Govan Manipulator 0.4mm × 0.2mm × 10.5mm Flattened Tip, 60° Angled. 75mm Handle</p> |  |

31 Choppers & Manipulators

| | | | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <p>CYSU1291 Spatula Ø0.5mm to Ø0.85mm × 10mm Spatula Tip, 45° Angled. 75mm Handle</p> |  | <p>CYSU1142 Button Manipulator 0.25mm × 0.3mm Button on 1mm Tip, 12mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1179 Mushroom Manipulator 0.5mm × 0.65mm Mushroom Tip, 10.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1428 Diablo Manipulator Ø0.35mm × 0.35mm Ball on 0.75mm Tip, 10.5mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1457 Barraquer Iris Spatula Ø0.65mm × 0.3mm Spatula Tip, 15mm 60° Angled Curved Shaft. 75mm Handle</p> |  | <p>CYSU1458 Lebuisson Manipulator 0.40mm × 0.50mm × 1.00mm Spatula Tip on Ø0.3mm Delicate Curved 12.6mm 90° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1432 Flat, Round Tip Manipulator Ø1.0mm × 0.15mm Round Tip, 10mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1434 Butcher Rotator Ø0.5mm Ball Tip, 10mm Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1208 Jaffe/Maltzman Manipulator Ø0.2mm × 0.7mm Tip Ø0.12mm × 0.65mm Extension at 50°, 12mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1207 Lester Hook Ø.4mm Ball on 0.8mm Finger, 11mm 45° Angled Shaft. 75mm Handle</p> |  |
| <p>CYSU1430 Osler Hook Ø0.2mm × 0.5mm Tip Ø0.15mm × 0.7mm Hook at 60° , 12.5mm 45° Angled Shaft. 75mm Handle</p> |  | <p>CYSU1109 'Y' Manipulator (Vertical) 2 × 0.55mm Tips with 0.8mm Spread, 45° Angled Shaft. 75mm Handle CYSU1204 - Horizontal</p> |  |









Choppers & Manipulators 32

Double-ended Choppers/Manipulators



| | | | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <p>CYSU1476 Mackool Iris Repositor 2 × 0.25mm Tips with 0.25mm Spread, 12.0mm 45° Angled Shaft. 75mm Handle</p> | | <p>CYSU1225 Kuglen Hook, Straight 0.25mm × 0.7mm Hook, 20mm Straight Shaft. 75mm Handle</p> | |
| <p>CYSU1433 Kuglen Hook, Angled 0.27mm × 1.1mm Hook, 10mm 45° Angled Shaft. 75mm Handle</p> | | <p>CYSU1112 Kuglen Hook 2 × Spiked Arms, 0.7mm Spread at Tips, 12mm 45° Angled Shaft. 75mm Handle</p> | |
| <p>CYSU1429 Trumpet Manipulator Ø0.7mm × 0.15mm Deep Trumpet Tip, 9.0mm 45° Angled Shaft. 75mm Handle</p> | | <p>CYSU1443 Micro Manipulator 0.6mm × 0.4mm Pointed Tip, Flattened 9.5mm 60° Angled Shaft. 75mm Handle</p> | |
| <p>CYSU1221 Phaco Chopper/Mushroom Manipulator 0.3mm × 1.5mm Chopper Tip 0.5mm × 0.65mm Mushroom Tip, 10.5mm 45° Angled Shafts. 52mm Handle</p> | | <p>CYSU1487 Phaco Chopper/Mushroom Manipulator, Fine 0.3mm × 0.95mm Chopper Tip 0.45mm × 0.65mm Mushroom Tip, 14.0mm/10.5mm 45° Angled Shafts. 52mm Handle</p> | |


33 Choppers & Manipulators

| | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>CYSU1238 Ursell Chopper/Manipulator 0.3mm × 1.5mm Tip, 12.5mm 45° Angled Shaft 2 × Ø0.2mm Tips at 120° Angle, 11.5mm 45° Angled Shaft. 52mm Handle</p>  | <p>CYSU1239 Diablo Chopper/Flat Button Manipulator Ø0.35mm × 0.35mm Ball on 0.75mm Tip Ø0.65mm × 0.35mm Button Tip, 10mm 45° Angled Shafts. 52mm Handle</p>  |
| <p>CYSU1431 Koch Chopper, 1.4mm/Button Manipulator 1.4mm Cylindrical Tip, 9.5mm 45° Shaft 0.25mm × 0.3mm Button Tip, 12mm 45° Angled Shaft. 52mm Handle</p>  | <p>CYSU1452 Ayliffe Chopper/Chopper 0.35mm × 1mm 45° Angled Blade, 10.5mm 45° Angled Shaft 0.28mm × 0.58mm Rectangular Prism Shaped Chopper, 10mm 40° Angled Shaft. 52mm Handle</p>  |
| <p>CYSU1248 Spatula/Horizontal 'Y' Manipulator Ø0.5mm to Ø0.85mm × 10mm Spatula Tip 2 × 0.55mm Tips with 0.78mm Spread, 9.5mm 45° Angled Shafts. 52mm Handle</p>  | <p>CYSU1496 Phaco Chopper/Govan Manipulator 0.3mm × 1.5mm Chopper Tip, Rounded Bottom, 11.5mm 45° Angled Shaft. 0.4mm × 0.2mm × 10.5mm Flattened Tip at 60° Angle. 70mm Handle</p>  |
| <p>CYSU1437 Koch Chopper/Manipulator (Rahman) 1.4mm Cylindrical Tip, 9.5mm 45° Angled Shaft. 0.2mm × 0.7mm × 13mm 60° Angled Shaft. 52mm Handle</p>  | <p>CYSU1438 Right Greens Chopper/Manipulator (Ikram) 1.4mm Chopper Tip, 10mm 45° Angled Shaft, 20° Bend. 0.2mm × 0.7mm × 13mm 60° Angled Shaft. 52mm Handle</p>  |

LASEK/PRK/Cross-Linking 34

CYSU1113

Corneal Epithelium Spatula, Flat
7.5mm × 0.2mm × 2.3mm Flat Blade.
75mm Handle



CYSU1187

Corneal Epithelium Spatula, Round
2.65mm × 5.5mm × 0.35mm Rounded Blade.
75mm Handle



CYSU1117

Epithelium Hoe, Flat
2.2mm × 0.45mm × 1.15mm Micro Hoe,
10.5mm 45° Angled Shaft.
75mm Handle



CYSU1185


Epithelium Hoe, Curved
1.7mm × 0.3mm × 1.0mm Micro Hoe,
10.5mm 45° Angled Shaft.
75mm Handle



Epithelium Wells

ID Well, Angled 45°. 75mm Handle


| | | |
|----------|-----------------|----------------|
| CYSU1201 | 6.5mm Shallow | 3.0mm × 6.5mm |
| CYSU1200 | 7.0mm Shallow | 3.0mm × 7.0mm |
| CYSU1199 | 7.5mm Shallow | 3.0mm × 7.5mm |
| CYSU1198 | 8.0mm Shallow | 3.0mm × 8.0mm |
| CYSU1197 | 8.5mm Shallow | 3.0mm × 8.5mm |
| CYSU1203 | 9.0mm Shallow | 3.0mm × 9.0mm |
| CYSU1132 | 9.5mm Shallow | 3.0mm × 9.5mm |
| CYSU1202 | 10.0mm Shallow | 3.0mm × 10.0mm |
| CYSU1196 | 9.0mm Trepining | 3.0mm × 9.0mm |



35 Miscellaneous



Vite Forcep 23g , 25g



BP HANDLE CYSU 300





Certificate of Registration

This certifies that the Quality Management System of

Huai'an Chaoyue Medical Supplies Factory

No.18-18, Industrial Park, No.18, Xi'an Road,
Huai'an City, Jiangsu Province, China 223021

has been assessed by NSF-ISR and found to be in conformance to the following standards(s):

ISO 13485:2016


Scope of Registration:

Manufacture of Non Sterile Non Active Ophthalmic Surgical Instruments

Exclusions: Product Design



| | |
|-------------------------|---------------|
| Certificate Number: | CNQMSMD038485 |
| Certificate Issue Date: | 02-JAN-2019 |
| Registration Date: | 02-JAN-2019 |
| Expiration Date*: | 01-JAN-2022 |



 Carl Blazik,
 Director- Technical Operations &
 Business Units, NSF-ISR, Ltd.

NSF International Strategic Registrations

789 North Dixboro Road, Ann Arbor, Michigan 48105 | (888) NSF-9000 | www.nsf-isr.org

Authorized Registration and /or Accreditation Marks. This certificate is property of NSF-ISR and must be returned upon request.
*Company is audited for conformance at regular intervals. To verify registrations call (888) NSF-9000 or visit our web site at www.nsf-isr.org

ISO CERTIFICATE





淮安超越医疗用品厂

CHAO YUE MEDICAL INSTRUMENTS FACTORY

T: +86 791 83814 302 / 517 83762 496

M: +86 189 7087 8605

Factory add: 18-18 , Xi'an Rd Industry Zone ,
Huai'an City, Jiangsu Province, PR.China

W: www.chaoyue-medical.com

E: sales@chaoyue-medical.com

